We can fill your needs.



We've
Expanded
Our Line Of
Central Line
And V.A.D. Access
Trays and Kits



Sterile Central Line Dressing Change Trays



Central Venous Catheter Dressing Change Tray with Tegaderm[™] Dressing and Isopropyl Alcohol Swabsticks

Provides the necessary components, including presaturated swabsticks and ointment, in sequence of use on a sterile field, for maintenance care of central venous catheters. The featured dressing is Tegaderm[™], a thin, flexible, transparent, polyurethane film coated with a hypoallergenic, acrylate adhesive. The defatting agent is 70% isopropyl alcohol on swabsticks, for those procedures not requiring the use of acetone.

Contents:

- · One trav
- · One CSR overwrap
- · One mask
- · One pair gloves
- One povidone iodine ointment
- · Three iodophor PVP swabsticks
- Three isopropyl alcohol swabsticks
- One 2"x 2" non-adhering dressing
- One 2" x 2" gauze dressing
- One 4" x 4" gauze dressing
- One 2" x 2" non-adhering split dressing
- One 4" x 5" Tegaderm™ transparent dressing
- · One roll medical tape
- · One dressing change label

CATALOG NO. 829

Packed 30 trays/carton

Central Venous Catheter Dressing Change Tray with Bioclusive® Dressing and Isopropyl Alcohol Swabsticks

Same as Cat. No. 829, with Bioclusive® transparent dressing in place of Tegaderm™ dressing.

CATALOG NO. 828

Packed 30 trays/carton



Central Venous Catheter Dressing Change Tray with Non-Adhering Dressings and Acetone Alcohol Swabsticks

Provides the necessary components, including presaturated swabsticks and ointment, in sequence of use on a sterile field, for maintenance care of central venous catheters. The featured dressings are highly absorbent and do not stick to wounds.

Contents:

- · One tray
- One CSR overwrap
- · One mask
- One pair gloves
- · One povidone iodine ointment
- · Two iodophor PVP swabsticks
- · Three acetone alcohol swabsticks
- · One benzoin swabstick
- One cotton tipped applicator
- One 2" x 2" non-adhering split dressing
- Three 2" x 2" non-adhering dressings
- One 3" x 4" non-adhering dressing (occlusive)
- · One roll medical tape
- · One dressing change label
- One polybag with twist tie

CATALOG NO. 825

Packed 30 travs/carton



Central Venous Catheter Dressing Change Tray with Bioclusive® Dressing and Acetone Alcohol Swabsticks

Provides the necessary components, including presaturated swabsticks and ointment, in sequence of use on a sterile field, for maintenance care of central venous catheters. The featured dressing is Bioclusive®, a thin, flexible, transparent, polyurethane film.

Contents:

- · One tray
- One CSR overwrap
- One mask
- One pair gloves
- One povidone iodine ointment
- · Three iodophor PVP swabsticks
- Three acetone alcohol swabsticks
- One cotton tipped applicator
- One 2" x 2" non-adhering dressing
- One 3" x 3" gauze dressing
- One 4" x 5" Bioclusive® transparent dressina
- One roll medical tape
- · One dressing change label
- · One polybag with twist tie

CATALOG NO. 827

Packed 30 trays/carton



Central Venous Catheter Dressing Change Tray with <u>Tegaderm</u>™ Dressing and Acetone Alcohol Swabsticks

Provides the necessary components, including presaturated swabsticks and ointment, in sequence of use on a sterile field, for maintenance care of central venous catheters. The featured dressing is Tegaderm™, a thin, flexible, transparent, polyurethane film coated with a hypoallergenic, acrylate adhesive.

Contents:

- · One tray
- One CSR overwrap
- One mask
- · One pair gloves
- · One povidone iodine ointment
- · Three acetone alcohol swabsticks
- · Three iodophor PVP swabsticks
- One 4" x 4" gauze dressing
- One 2" x 2" non-woven split dressing
- One 2" x 2" non-adhering dressing
- One 2" x 2" gauze dressing
- One 4" x 5" Tegaderm™ transparent dressing
- One roll medical tape
- · One dressing change label

CATALOG NO. 833

Packed 30 trays/carton





Central Venous Catheter Dressing Change Kit with <u>Tegaderm</u>™ Dressing and Skin Protectant Swabsticks

Provides the necessary components to streamline maintenance care of the central venous catheter. Components are assembled in sequence of use in a flexible package and include presaturated swabsticks and ointment. The featured dressing is Tegaderm™, a thin, flexible, transparent polyurethane film dressing. Instead of a tray, this kit comes soft packed with a polylined sterile field.

Contents:

- · One polylined drape
- · One pair gloves
- · One iodophor PVP ointment
- Three isopropyl alcohol swabsticks
- · Three iodophor PVP swabsticks
- One skin protectant swabstick
- One 2" x 2" non-woven dressing
- One 2" x 2" non-woven split dressing
- One 3" x 3" gauze dressing
- One 4 " x 5 " Tegaderm™ transparent dressing
- · One roll medical tape
- · One dressing change label

CATALOG NO. 832

Packed 50 trays/carton





Central Line Dressing Change Kit with <u>Tegaderm</u>™ Dressing and Skin Protectant Swabsticks

Provides the necessary components to streamline maintenance care of the central line catheter. Components are assembled in sequence of use in a flexible package and include presaturated swabsticks and ointment. The featured dressing is Tegaderm™, a thin, flexible, transparent polyurethane film dressing. Instead of a tray, this kit comes soft packed with a polylined sterile field.

Contents:

- · One polylined drape
- · One mask
- One pair gloves
- Three isopropyl alcohol swabsticks
- · Three iodophor PVP swabsticks
- · One skin protectant swabstick
- One 2" x 2" non-woven dressing
- One 2" x 2" gauze dressing
- One 3" x 3" gauze dressing
- One 4" x 5" Tegaderm™ transparent dressing
- · One roll medical tape

CATALOG NO. 834

Packed 50 trays/carton





V.A.D. Access Tray with Tegaderm™ Dressing, Heparin and Saline

This kit can be used for the care of all types of vascular access devices, especially implantable access devices. It can be used prior to and following infusion treatments such as TPN, chemotherapy, antibiotic therapy and fluid replacement, as well as for maintenance care and frequent blood drawing.

Contents:

- One tray
- · One earloop mask
- One CSR overwrap
- · One pair gloves
- · Three isopropyl alcohol swabsticks
- · Three iodophor PVP swabsticks
- · One skin protectant swabstick
- One 3"x 3" gauze dressing
- One 2" x 2" gauze dressing
- One 4" x 5" Tegaderm™ transparent dressing
- One roll medical tape
- · One injection cap
- Two 20g x 1" needles
- · One 20cc luer lock syringe
- · One 10cc luer lock syringe
- One 30cc sodium chloride 0.9%
- One 10cc heparin, 100 u/cc

CATALOG NO. 831

Packed 20 trays/carton

CUSTOM KITS



You pick the tray, we'll make it.

Don't see the tray you're looking for? Then why not custom design your own?

If our standard trays cannot be adapted to your particular procedure, we'd be pleased to custom design and manufacture one to your exact specifications, with the same high quality components found in our standard kits and trays.

With a custom tray, you can also specify the sequential packing of the components for proper aseptic technique and maximum efficiency. And since all components necessary to perform the procedure are prepared for you in a single, pre-sealed, sterile disposable tray, you can cut down on time spent ordering supplies and preparing for the procedure.

Sterile, disposable procedure trays are economical alternatives, saving your facility both time and money...not to mention the added benefit of simplified patient charges.

If you have an interest in exploring custom kit options, ask your Busse representative, or call our Customer Service Department at (800) 645-6526. Minimum quantities apply.

In this rapidly evolving environment of health care delivery, Busse is your prime advantage.

