### **Frazier Suction Tips**

# Get a Better Grip.

At Busse, we know how important the instrument is in a delicate procedure. We have redesigned our frazier suction tips to meet your surgical needs and exceed your expectations.

CATALOG NO. 162 - 8 FR

**CATALOG NO. 163 - 10 FR** 

**CATALOG NO. 164 – 12 FR** 

Packaged Sterile
5 dispenser boxes of 10 each

#### Sterile obturator

clears obstructions and allows the malleable shaft to be bent to the desired position without occlusion.

#### Longer tube

for easier access to hard to reach areas.

#### Anodized aluminum shaft

provides flexibility during procedures. Available in three sizes: 8FR, 10FR and 12FR for procedure versatility.

#### Control vent

is conveniently located for continuous or intermittent suctioning.

#### Ergonomically designed handle

is longer, allowing for better dexterity. The square shape ensures non-slip grip.

#### Universal connector

facilitates easy connection of tubing.

## DUSSE Hospital Disposables

P.O. Box 011067 Hauppauge, NY 11788-0920 800.645.6526 • Fax: 631.435.4721 www.busseinc.com • bussemail@aol.com

© 2000 Busse Hospital Disposables E.O.E. - M/F/D/V