Personal Protection Gown

Protection, Comfort and Convenience for the Caregiver.

These gowns are the ideal impervious barrier protection for caregivers in the emergency room as well as in patient rooms. Thumbhooks prevent the gown from riding up and exposing wrists. The pullover design eliminates ties at the neck and the perforations in the back allow for fast removal. Gowns come in universal size and are made of latex-free, blue embossed polyethylene.



CATALOG NO. 235 Individually packed 15/box, 5 boxes/carton



THUMBHOOK STIRRUPS
Makes glove donning easier and protects against wrist exposure.



COMFORTABLE & CONVENIENT Open back keeps wearer cool and comfortable. Perforations allow for fast, easy removal.



Individually packed 15/box.

Hooks on to mounting bracket on most personal protection stations.

