

NEW ADDITION TO OUR PAIN MANAGEMENT LINE!



PUMP REFILL KITS FROM BUSSE

Convenient • Cost Effective • Universal

Busse's Pump Refill Kit:

- Provides quality components needed in one convenient pack
- Is specifically designed to be universal (with user supplied template)
- Simplifies inventory management and optimizes purchasing
- Offers drug-free packaging for more flexible storage
- Is available through all distribution channels



Busse Pump Refill Kit

Cat No. 668

Sterile Packaged 10 per case

Components:

- 1 ea. Tray
- 1 ea. Fenestrated drape
- 1 pr. Walleled exam gloves
- 1 ea. 20mL syringe
- 2 ea. Gauze sponges
- 1 ea. Huber needle 22g x 1.5"
- 1 ea. Huber needle 22g x 2"
- 1 ea. Extension set
- 1 ea. Micron filter
- 1 ea. Iodophor PVP swabsticks (3/pk)
- 2 ea. Alcohol prep pad
- 2 ea. Adhesive bandage

Customized to your needs

If our standard tray does not meet all of your exact requirements, we can customize one for you. Please call our Customer Service Department at **1-800-645-6526** for our sales representative in your area.

NOT MADE WITH
Natural Rubber Latex



Please go to our website or scan the QR code on your mobile device to see our full line of kits and trays.

STERILE EO



ISO 13485:2003
Certified
FM 77400

busse
Hospital Disposables

75 Arkay Drive, Hauppauge, NY 11788-0920
Tel: 800.645.6526 Fax: 631.435.4721
www.busseinc.com mail@busseinc.com
©2016 Busse Hospital Disposables
Product specifications are subject to change.
EEO / AA