SMS Protective Apparel

TIMES BETTER

When protecting your staff is as critical as protecting your patients, choose SMS apparel from Busse. You'll be providing more resistance to blood strike through as well as comfort with this extra strong material.

BETTER PROTECTION

More fluid resistant than regular spunbonded material, the multi-layered SMS apparel offers better protection for any level of exposure in a variety of health care situations.

BETTER COMFORT

Sized generously and softer than regular spunbonded materials, latex free SMS apparel drapes easily, making gowns, labcoats and coveralls more comfortable.

BETTER VALUE

Once again, Busse gives you the quality and technology that others offer, but at a fraction of the cost.



SMS is made of meltblown polypropylene sandwiched between two layers of spunbonded polypropylene to create a stronger, yet more comfortable barrier.

LATEX FREE



ISOLATION GOWNS

Exposure Anticipated: potential for light fluid contact and contact isolation.

- Elastic cuffs
- Tape back neck closure
- Full back
- Waist ties
- Gown is sewn, not glued
- · One size fits all
- Yellow

CATALOG NO. 242

Packed 100/case (10/bag, 10 bags/case)

COVERALLS

Exposure Anticipated: potential for moderate fluid contact.

- · Open cuffs
- Blue

CATALOG NO. 245 — Large CATALOG NO. 246 — X-Large

Packed 24/case

- · Elastic wrist cuffs
- Open ankle cuffs
- Back elastic gather
- White

CATALOG NO. 247 — Large CATALOG NO. 248 — X-Large Packed 24/case

LAB COATS

Exposure Anticipated: potential for moderate fluid contact and extended use.

- Knit cuffs
- Traditional collar
- 6 snap front closures
- 3 reinforced pockets
- Back vent
- Neck loop for hanging
- White

CATALOG NO. 222 – Small/Medium CATALOG NO. 225 – Large/X-Large

CATALOG NO. 229 - XX Large

Packed 25/case

© 2002 Busse Hospital Disposables E.O.E. - M/F/D/V

P.O. Box 11067 Hauppauge, NY 11788-0920 800.645.6526 • Fax: 631.435.4721 • www.busseinc.com • mail@busseinc.com