Umbilical Cord Clamps and OB Kits



Busse has the answer for every delivery situation.

Busse offers easy-to-use Posi-Grip™ Umbilical Cord Clamps for use in the delivery room and complete obstetrical kits for use by paramedics and police out in the field.

See reverse for more information.

BUSSE DELIVERS QUALITY





Umbilical Cord Clamps and OB Kits



POSI-GRIP™ UMBILICAL CORD CLAMP

ONCE IT'S CLOSED, IT'S LOCKED!

The single-use clamp is designed to seal the umbilical cord, preventing loss of blood as the cord dries beneath it. It can be easily placed with one hand, even if gloved and slippery. The interlocking teeth hold the cord securely and keep it stationary. Once the protected closure is in place and closed, it cannot open accidently. Works with any style cord cutter.

CATALOG NO. 384



Sterile - Individual Packets, 128/carton (2 boxes of 64 each)

CATALOG NO. 385



Non-sterile - Bulk Packed 250/carton

OBSTETRICAL KITS

A complete disposable obstetrical kit available in cardboard box or soft pouch packaging. Key components are individually packed sterile for emergency delivery situations.

Items included in both pouch and boxed versions:

- · One pair sterile exam gloves
- · One disposable plastic apron
- One 17" x 24" underpad
- One 23" x 36" receiving blanket
- One 48" x 36" drape sheet
- Four absorbent towels
- Two OB towelettes
- Four sterile 4" x 3" gauze sponges
- One sterile #10 scalpel
- · One sterile OB Pad
- One sterile 2oz bulb syringe
- · Two alcohol prep pads
- Two sterile Posi-Grip™ umbilical cord clamps

P.O. Box 11067 Hauppauge, NY 11788-0920

- · One clear polybag with tie
- · One white polybag with tie



POUCH NO. 388

Packed 10 kits/carton



Packed 10 kits/carton

